


Perfect™ Tricolor Protein Ladder

Cat. No.	size
E3210-01	50 loadings

Storage Conditions: Store at -20°C.




Prestained Protein Ladder, composed of precisely spaced recombinant bands, for sizing polypeptides ranging from 10 kDa to 245 kDa.

Description:

- The Perfect™ Tricolor Protein Ladder has twelve bands corresponding to protein size: 11, 17, 20, 25, 35, 48, 63, 75, 100, 135, 180 and 245 kDa.
- The 25 kDa band is covalently coupled to a green dye while the 75 kDa band is covalently coupled to an orange dye. The remaining bands are covalently coupled to a blue dye.
- Comes in a ready-to-use format.

Loading:

The recommended amount of size marker to load on a mini-gel is 5 µl per lane.



Tris-Glycine Gradient gel 4-20%

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23